

CR PLATES AND ACCESSORIES



AIR BULB BRUSH

Air Bulb Brush.

- *Ideal for removing dust on object in macro photography.*
- *High-quality rubber.*
- *Great way to remove dust off sensitive camera parts without physical contact.*

[Read More](#)

CODE SC300AIR

Categories: [Cleaning Products](#), [CR Plates and Accessories](#), [Infection Control](#)

Tag: [Cleaning accessories](#)



ATX CR PLATE CLEANER KIT

ATX CR Plate Cleaner Kit:

- *1 x 250ml bottle ATX CR Screen Cleaner*
- *1 x Air Bulb Brush*
- *1 x ATX Equipment Wipes (75 pack)*

[Read More](#)

CODE SC300KIT

Category: [CR Plates and Accessories](#)

Tags: [Cleaning Solution](#), [CR plate cleaner](#)



X-RAY PLATE DISPOSABLE SMALL PLASTIC COVERS (PACK OF 25)

X-ray Plate Protector. Small Plastic Covers, sold per pack of 25 bags.

[Read More](#)

CODE DZ3336

Categories: [Cleaning Products](#), [CR Plates and Accessories](#)

Tags: [disposable](#), [Plastic covers](#), [protection](#), [X-ray plate protection cover](#)



X-RAY PLATE DISPOSABLE LARGE PLASTIC COVERS (PACK OF 25)

X-ray Plate Protector. Large Plastic Covers, sold per pack of 25 bags.

[Read More](#)

CODE DZ5050

Categories: [Cleaning Products](#), [CR Plates and Accessories](#)

Tags: [disposable](#), [Plastic covers](#), [protection](#), [X-ray plate protection cover](#)



AGFA PLATE CLEANER (500MLS PER BOTTLE)

AGFA Plate Cleaner (500ml per bottle)

CR Phosphor Plate Cleaning is an easy-to-use, ethanol-based cleaner that ensures longer life and higher performance for CR Imaging Plates.

[Read More](#)

CODE SC300

Categories: [Cleaning Products](#), [CR Plates and Accessories](#)

Tag: [Cleaning Solution](#)



ATX CR PLATE CLEANER

ATX CR Plate Cleaner, is easy to use; Plate Cleaning Solution 250mls per bottle.

CR X-ray Plates are highly sensitive and requires regular cleaning. We recommended cleaning your plate every 200 scans, or every 30 days (depending of your work environment).

[Read More](#)

CODE SC300ATX

Categories: [Cleaning Products](#), [CR Plates and Accessories](#), [Infection Control](#)

Tags: [Cleaning Solution](#), [CR plate cleaner](#)