


DENTAL IMAGING SYSTEMS



RAPID PARALLELING KIT WITH BITE WING

RAPiD is a right angle positioning device that assures correct film/cone alignment and prevents "cone cutting". It can be used for anterior, posterior, and bite-wing x-ray exams using any brand and size (3, 2, 1, and 0) of dental film.

[Read More](#)

CODE 408560

Categories: [Dental Imaging Systems](#), [Positioning Tools](#)

Tags: [dental](#), [positioning kit](#), [Smart](#)



RAPID ENDODONTIC POSITIONER

Includes 1 Aiming Ring, 1 Indicator Arm, and 3 Bite Blocks.

[Read More](#)

CODE 40303

Categories: [Dental Imaging Systems](#), [Positioning Tools](#)

Tags: [endodontic](#), [intraoral X-Ray](#), [positioning kit](#)



ENDOGRAPH DC MOBILE VERSION

Endograph DC matches any practice, as it is available with 3 different sizes of extension arms that can be adapted to virtually any room. Upon request, Endograph DC can come also as a mobile unit mounted on trolley, particularly indicated in emergency contexts. Additionally, the standard availability of preset settings for films, Photostimulable Phosphor Plates (PSP), or Digital Detectors, allows the use with any receptor, from any manufacturer.

[Read More](#)

CODE 9459000313

Categories: [Dental Equipment](#), [Dental Imaging Systems](#), [Intraoral X-ray](#)

Tags: [intraoral X-Ray](#), [Mobile](#), [wireless](#)



ENDOGRAPH DC INTRAORAL X-RAY

- **The innovation of a breakthrough user interface**

Endograph DC stands out completely from conventional intraoral units and makes a step further in the direction of user friendliness: its beautiful, shield-shaped control unit integrates a full-size LCD display showing the parameters efficiently organized in working areas and, for the first time, cleverly combining icons and tooth name. Additionally, the user interface can be customized in many languages, thus helping the operator to establish a communication channel with the patients. The levels of brightness and contrast can be

Tailored veterinary system and solutions including digital - CR - DR x-ray, other diagnostic equipment supply and service.

adjusted, depending on the environment where the unit is being installed.

- **Wireless X-ray pushbutton**

Endograph DC provides wireless control exposure. The remote timer can be mounted as well out of the way of the beam, to combine settings and exposure out of the patient area.

Huge range of configurations

Endograph DC matches any practice, as it is available with 3 different sizes of extension arms that can be adapted to virtually any room. Upon request, Endograph DC can come also as a mobile unit mounted on trolley, particularly indicated in emergency contexts. Additionally, the standard availability of preset settings for films, Photostimulable Phosphor Plates (PSP), or Digital Detectors, allows the use with any receptor, from any manufacturer.

[Read More](#)

CODE N/A

Categories: [Dental Equipment](#), [Dental Imaging Systems](#), [Intraoral X-ray](#)