

DIGITAL RADIOGRAPHY



LONGTAIL DR SYSTEM - 108X43CM

High performance cassette sized DR panel

The worlds only large single plate DR system.

- Huge 108 x 43cm field for large image acquisition.
- High DQE - to achieve outstanding image quality at low dose
- Ultra-sensitive AED - eliminates generator integration
- Lowest dose
- Fast - full resolution image viewable in a few seconds
- Low maintenance - frequent re-calibration not required
- Battery for wireless option with 10 hours of standby time
- No external interface or power box

for more information please visit www.longtaildr.com

[Read More](#)

CODE 2ATX-DRGadox W 14x17-1

Category: [Digital Radiography](#)



ATX DR SYSTEM - 35X43CM

High performance cassette sized DR panel

- High DQE - to achieve outstanding image quality at low dose
- Ultra-sensitive AED - eliminates generator integration
- Lowest dose
- Fast - full resolution image viewable in a few seconds
- Low maintenance - frequent re-calibration not required
- Battery for wireless option with 10 hours of standby time
- No external interface or power box

[Read More](#)

CODE 2ATX-DRGadox W 14x17

Category: [Digital Radiography](#)



KONICA MINOLTA CR - SIGMA II

COMPLETE SYSTEM WITH ALL-IN-ONE TOUCH SCREEN COMPUTER AND IMAGE PILOT SOFTWARE.

The REGIUS Σ is Konica Minolta's smallest and lightest CR. The desktop CR has a footprint of only 0.31m² and weighs a mere 28kg.

ENVIRONMENTALLY FRIENDLY



The REGIUS Σ was designed with the environment in mind. Energy consumption has been reduced to 100VA - about the same as a household light bulb. That's one tenth the power of the smallest film processor. Important with power costs soaring.

HIGH RESOLUTION



The REGIUS Σ supports pixel sizes of 87.5 μ m and 175 μ m for all cassette sizes, achieving a maximum resolution of 4020x4892 pixels (35x43cm at 87.5 μ m).

INTEGRAL PROCESSING

Integral Processing - the result of Konica Minolta's vast imaging experience - is a robust algorithm that produces the optimum diagnostic image automatically. This patented technique greatly simplifies x-ray image acquisition while delivering incredibly consistent image quality and virtually eliminating retakes.



[Read More](#)

CODE 2KMRegius CR

Category: [Digital Radiography](#)

Tailored veterinary system and solutions including digital - CR - DR x-ray, other diagnostic equipment supply and service.



ATX DR SYSTEM - 43X43CM

High performance cassette sized DR panel

Now with new improved acquisition software for optimal image processing.

- **High DQE - to achieve outstanding image quality at low dose**
- **Ultra-sensitive AED - eliminates generator integration**
- **Lowest dose**
- **Fast - full resolution image viewable in a few seconds**
- **Low maintenance - frequent re-calibration not required**
- **Battery for wireless option with 10 hours of standby time**
- **No external interface or power box**

[Read More](#)

CODE 2ATX-DRGadox W 17x17

Category: [Digital Radiography](#)