

# SHEETS AND GOWNS

## DISPOSABLE FITTED WHITE SHEET BOX 50

Disposable Fitted White Sheet. Elastic all around (Anti-static material)  
Cartons of 50 sheets.

Our range of disposable bed sheets are strong breathable, with an anti-static soft construction to maximise comfort for your patients.

[Read More](#)

**CODE** C-DIS1003

**Category:** [Sheets and Gowns](#)

**Tag:** [Disposable Fitted White Sheet](#)



## KN95 RESPIRATOR



CE Certified, KN95 (P2, P95 equivalent) - not TGA Registered.

Exceeds minimum filtering standards for micro particles for optimum protection.

Mask material is able to filter particles as small as 300nm or 0.3 micron.

This link is a guide to various particle sizes.

[Visualise Relative Size of Particles](#)

**Priced from : \$4.20 Each**

**QTY: 5, 10, 50**

Individually wrapped

**Comfortable**

Elastic ear-loops with no pressure to the ears.

**Skin-friendly**

Made of non-woven material, they are featuring good skin tolerance and comfortable to wear. Without glass fibers.

**Breathable**

Safe, breathable, and hypoallergenic, high filtration capacity.


One size fits most people. Great choices for you or your family.

[Read More](#)

**CODE** C-DIS1012.

**Categories:** [Gloves](#), [Infection Control](#), [Sheets and Gowns](#)

**Tags:** [face mask](#), [KN95](#), [Medical](#), [P2](#), [PN95](#), [Respirator](#), [ZHIYI](#)





## PAPER BED ROLL 2PLY

*2Ply perforated quality paper with a universal fit to protect your treatment table. **With 3 Different Sizes of Paper Roll.***

- *Small*
- *Medium*
- *Large*

[Read More](#)

**CODE** C-DIS1015

**Category:** [Sheets and Gowns](#)

**Tag:** [Paper Roll](#)

## DISPOSABLE PLASTIC PILLOW COVER



Disposable Plastic Pillow Cover. Strong light blue vinyl (wipe clean and washable) sold individually.

[Read More](#)

**CODE** C-DIS1008

**Category:** [Sheets and Gowns](#)

**Tag:** [Disposable Plastic Pillow Cover](#)



## DISPOSABLE WHITE PILLOW COVER

*Disposable White Pillow Cover, strong, anti-static material, but flexible and easy to remove and replace. Carton of 100 pillows.*

Size: 48cm x 70cm

[Read More](#)

**CODE** C-DIS1007

**Category:** [Sheets and Gowns](#)

**Tags:** [Covers](#), [disposable](#), [Disposable White Pillow Cover](#), [Pillow](#), [White](#)





## DISPOSABLE FITTED SHEET BOX

Disposable Fitted White Sheet. Elastic all around (Anti-static material)  
Cartons of 70 sheets. Our range of disposable bed sheets are strong  
breathable, with an anti-static soft construction to maximize comfort for your  
patients. Flat Sheets are now available in both blue and white colors. Also  
available in light blue.


Size: 65cm x 190cm

[Read More](#)

**CODE** C-DIS1004

**Category:** [Sheets and Gowns](#)

**Tags:** [disposable](#), [Disposable Fitted White Sheet](#), [Flat](#), [Sheets](#), [White](#)



## EXTRA LARGE PATIENT GOWNS (100)

A simple change can have a large impact on patient dignity.

- Reliable Quality 55gsm Polyethylene breathable and anti-static material
- Sleeveless blue gown
- Sewn-on wrap-around tie for patient comfort and gown resilience
- Non transparent for patient confidence, privacy and dignity
- Over locked edges for durability and gown strength
- Cost effective

[Read More](#)

**CODE** C-DIS1002

**Category:** [Sheets and Gowns](#)

**Tag:** [Disposable XL Patient Gown](#)



## DISPOSABLE HEADREST COVER (BOX OF 1000)

*Disposable Hygienic Headrest Cover with cross cut (45gsm) - available in boxes of 1000.*

### **FEATURES & BENEFITS**


- *The Covers are clean, hypoallergenic and provide a reliable sanitary environment for your clients.*
- *Strong, breathable and anti-static with soft construction make them super comfortable for your clients.*
- *100% recyclable for the environmentally conscious.*
- *Quick and convenient for busy therapist, making bacteria control a breeze.*

[Read More](#)

**CODE** C-DIS1011

**Category:** [Sheets and Gowns](#)

**Tag:** [face-rest covers](#)



## DISPOSABLE HEADREST COVER - CIRCLE CUT

*Disposable Hygienic Headrest Cover circle cut (35gsm) - available in boxes of 1000.*

### ***FEATURES & BENEFITS***

- *The Covers are clean, hypoallergenic and provide a reliable sanitary environment for your clients.*
- *Strong, breathable and anti-static with soft construction make them super comfortable for your clients.*
- *100% recyclable for the environmentally conscious.*
- *Quick and convenient for busy therapist, making bacteria control a breeze.*

[Read More](#)

**CODE** C-DIS1010

**Category:** [Sheets and Gowns](#)

**Tag:** [face-rest covers](#)



## DISPOSABLE PLASTIC PILLOW COVER WITH ZIP

Disposable Plastic Pillow Cover with Zip. strong light blue vinyl (wipe clean and washable) sold individually.

[Read More](#)

**CODE** C-DIS1009

**Category:** [Sheets and Gowns](#)

**Tag:** [Disposable Plastic Pillow Cover with Zip](#)



## DISPOSABLE DARK BLUE PILLOW COVERS

Disposable Blue Pillow Covers

We know it is important to provide clients with a comfortable environment while keeping your equipment clean and sanitary.

Our Pillow covers are strong and anti-static material. Carton of 100 pillows.

[Read More](#)

**CODE** C-DIS1006

**Category:** [Sheets and Gowns](#)

**Tag:** [Disposable Blue Pillow Cover](#)



## DISPOSABLE FLAT WHITE SHEET

Disposable Flat White Sheet,

2400x700mm and made with 45gsm,

(anti-static material) Cartons of 100 sheets.

[Read More](#)

**CODE** C-DIS1005

**Category:** [Sheets and Gowns](#)

**Tags:** [Bed Sheet](#), [disposable](#), [Disposable Flat White Sheet](#), [Disposable Linen](#), [flat sheet](#), [Infection Control](#), [white sheet](#)

## DISPOSABLE PATIENT GOWN (BOX OF 100)

A simple change can have a large impact on patient dignity.

- Reliable Quality 55gsm Polyethylene breathable and anti-static material
- Sewn-on wrap-around tie for patient comfort and gown resilience
- Non transparent for patient confidence, privacy and dignity
- Over locked edges for durability and gown strength
- Cost effective
- (colour may vary from previous stock)

[Read More](#)

**CODE** C-DIS1001

**Category:** [Sheets and Gowns](#)

**Tags:** [Blue](#), [disposable](#), [Disposable patient gowns](#), [Gowns](#), [Patient](#)

