

# SCALING TOOLS

## BARIUM PASTE 20ML

Barium Sulphate, radiopaque paste that provides excellent radiographic visualization of the true hoof wall. The tube contains 20mls of paste. The paste appears as an opaque (white line) on the image, when applied in a continuous bead from the coronary band to the ground surface. Due to the very thin depth of horn at the most proximal crown of the hoof wall the actual horn margin is rarely visible on even the softest penetrating views.

[Read More](#)

**CODE** C-GEN1010

**Categories:** [Contrasts](#), [Lead Letters](#), [Markers](#), [Scaling Tools](#), [Skin Markers](#)

**Tags:** [barium](#), [contrast](#), [equine](#), [hoof](#), [laminitis](#), [lateromedial](#)





## AKUCAL SCALING TOOL ( 30MM REPLACEMENT MARKER )

30mm

Replacement stainless steel marker for Image Scaling Tool

EACH

[Read More](#)

**CODE** E-ACC1032-30

**Category:** [Scaling Tools](#)

**Tags:** [Akucal](#), [Calibration Tools](#), [flexible](#), [image scaling device](#), [positioning guide](#), [replacement marker](#), [scaling](#), [Scaling tool](#), [surgical grade](#)



## AKUCAL SCALING TOOL ( 25MM REPLACEMENT MARKER )

25mm

Replacement stainless steel marker for Image Scaling Tool

EACH

[Read More](#)

**CODE** E-ACC1030-25

**Category:** [Scaling Tools](#)

**Tags:** [Akucal](#), [Calibration Tools](#), [flexible](#), [image scaling device](#), [positioning guide](#), [replacement marker](#), [scaling](#), [Scaling tool](#), [surgical grade](#)



## AKUCAL IMAGE SCALING DEVICE ( 30MM ) (COMPLETE)

25mm image scaling device

COMPLETE with all parts including-

- Suction cup base for Horizontal or Vertical use
- Flexible arm to position marker at anatomy level
- Aluminium tip avoids breakage from wear and tear
- Numerical size indicators visible in every image ( 25mm, 25.4mm, 30mm )
- Surgical Grade stainless steel marker
- Lever to engage or disengage suction

[Read More](#)

**CODE** E-ACC1031-30

**Category:** [Scaling Tools](#)

**Tag:** [Akucal](#)



## AKUCAL IMAGE SCALING DEVICE ( 25MM ) (COMPLETE)

25mm image scaling device

COMPLETE with all parts including-

- Suction cup base for Horizontal or Vertical use
- Flexible arm to position marker at anatomy level
- Aluminium tip avoids breakage from wear and tear
- Numerical size indicators visible in every image ( 25mm, 25.4mm, 30mm )
- Surgical Grade stainless steel marker
- Lever to engage or disengage suction

[Read More](#)

**CODE** E-ACC1029-25

**Category:** [Scaling Tools](#)

**Tag:** [Calibration Tools](#)

