

# RADIATION PROTECTION

## RADIATION PROTECTION BARRIER

Mobile or fixed lead barriers with Lead glass window - manufactured to size.

[Read More](#)

**CODE** N/A

**Categories:** [Lead Glass and Mobile Shields](#), [Personal Radiation Protection](#)






## CLOSED LOOP APRON HANGER

[Read More](#)

**CODE** C-PRO1014

**Categories:** [Personal Radiation Protection](#), [Radiation Protection](#)

**Tags:** [Apron Hanger](#), [Hanger](#), [Radiation Protection](#), [Steel Hanger](#)



## HEAVY-DUTY APRON HANGER

[Read More](#)

**CODE** C-PRO1013

**Categories:** [Personal Radiation Protection](#), [Radiation Protection](#)

**Tags:** [Hanger](#), [Lead Gown](#), [Radiation Protection](#), [Steel Hanger](#)



## WALL MOUNT SWING ARM APRON RACK

Wall Mount Swing Arm Apron Rack features 5 individual chrome shoulder shape swing arms that fold flush against the wall. Available in right and left folding models. Simply hang two models side-by-side for a 10 hanger wall console unit.

[Read More](#)

**CODE** C-PRO1005

**Category:** [Personal Radiation Protection](#)

**Tags:** [Gown](#), [Lead Gown](#), [Lead gown rack](#), [Rack](#)

## LEAD THYROID COLLAR

**Quick Release Collar** offers the same look and comfort as our model 621U, but without the hook & loop closure. Rear buckle closure. Available in our wide selection of nylons, ripstops and designer fabrics.

- **Standard Protection:** 0.5mm Pb Equivalency

**New Style Thyroid Collar** offers a trim yet shielded and comfortable design. Includes binding and rear hook & loop closure.

- **Standard Protection:** 0.5mm Pb Equivalency

Please note: Colour is light blue

[Read More](#)

**CODE** C-PRO1003-L

**Category:** [Personal Radiation Protection](#)

**Tags:** [lead protection](#), [Thyroid Collar](#)





## LEAD GOWN

Choose your Size

(XS) X-Small: 18" / 46 cm Wide - Fits Chest 30-34" / 76.2-86.4cm

(S) Small: 20" / 51 cm Wide - Fits Chest 34-38" / 86.4-96.5cm

(M) Medium: 22" / 56 cm Wide - Fits chest 38-42" / 96.5-106.7cm

(L) Large: 24" / 61 cm Wide - Fits Chest 42-46" 106.7-116.8cm

Regular lead = 0.5mm

[Read More](#)

**CODE** C-PRO1000

**Category:** [Personal Radiation Protection](#)

**Tags:** [Lead Gown](#), [lead shielding](#), [xray](#), [xray gown](#)



## NOLAN FILTER, GONAD SHIELD SYSTEM SET

Set of male and female gonad shields  
and belt.

[Read More](#)

**CODE** NFAPSET

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## NOLAN FILTER, GONAD SYSTEM BELT

Velcro belt to suspend gonad shields

[Read More](#)

**CODE** NFAPSBELT

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)





## NOLAN FILTER, A-P FEMALE GONAD SHIELD

A-P Female Gonad Shield for AP shots.

[Read More](#)

**CODE** NFAPFGS

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## NOLAN FILTER, A-P MALE GONAD

A-P Male Gonad Shield for AP shots.

[Read More](#)

**CODE** NFAPMGS

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## NOLAN FILTER, LATERAL SPINOUS

Lateral Spinous: Reduces scatter and skin dose and  
balances plate dose by blocking x-rays to the spinous process.

[Read More](#)

**CODE** NFLS

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## **NOLAN FILTER, A-P PARASPINAL BLOCKER (TOTAL)**

A-P Paraspinal Blocker (Total): Totally blocks anatomy either side of the spine.

[Read More](#)

**CODE** NFAPSBF

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## NOLAN FILTER, A-P PARASPINAL BLOCKER (PARTIAL)

A-P Paraspinal Blocker (Partial): Partially blocks anatomy either side of the spine.

[Read More](#)

**CODE** NFAPSBP

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## NOLAN FILTER, LATERAL GONAD

Lateral M/F Gonad shield for lateral exposures

[Read More](#)

**CODE** NFLG

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## NOLAN FILTER, LATERAL LUMBAR LUNG

Lateral Lumbar Lung: Lumbar Lung shield for lateral exposures

[Read More](#)

**CODE** NFLLL

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## LATERAL STERNUM

Lateral Sternum: Shield for lateral exposures

[Read More](#)

**CODE** NFSB

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)





## NOLAN FILTER, LATERAL THORACIC

Lateral Thoracic. Compensating wedge

[Read More](#)

**CODE** NFLT

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## NOLAN FILTER LATERAL CERVICO THORACIC


Accommodates relative thickness of the shoulder to the Cervical or Thyroid region.

[Read More](#)

**CODE** NFLCT

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## NOLAN FILTER FACE PLATE

Mounts directly on to the collimator to  
accommodate multiple filters

[Read More](#)

**CODE** NFFP

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)



## NOLAN FILTERS, A-P CERVICO THORACIC

A-P Cervico Thoracic Compensating wedge

[Read More](#)

**CODE** NFAPCT

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#)



## NOLAN FILTER POINT FILTERS

***POINT FILTERS. Can be used for compensation in all views.***


Adjust the exposure to the place according to anatomical thickness.  
Lower until the bottom is just above the iliac crests.

[Read More](#)

**CODE** NFPP

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)



## OCTOSTOP HALF INGOT FILTER

Half Ingot Filter


L29 x W19 x H4.5 cm approx 1.8kg.

[Read More](#)

**CODE** C-FIL1034

**Category:** [Collimator Filters](#)

**Tags:** [boomerang](#), [compensating filter](#), [filter](#), [Filters and Grids](#), [ingot](#)



## MOBILE APRON RACK 1.17MT

[Read More](#)

**CODE** C-PRO1009

**Category:** [Personal Radiation Protection](#)

**Tags:** [Hangers](#), [Racks & Shielding](#)



## LEAD GLOVES

[Read More](#)

**CODE** PSG35

**Category:** [Radiation Protection](#)

**Tag:** [Uncategorised](#)





## VET MITT

[Read More](#)

**CODE** C-PRO1002

**Category:** [Personal Radiation Protection](#)

**Tag:** [Lead Gloves](#)



## OCTOSTOP CHAINLESS BOOMERANG

Regular Boomerang with no chain. **PRICE: TBA**

[Read More](#)

**CODE** OCT01NC

**Category:** [Radiation Protection](#)

**Tags:** [boomerang](#), [compensating filter](#), [filter](#), [Filters and Grids](#)



## OPEN PALM SHIELD

[Read More](#)

**CODE** C-PRO1001

**Category:** [Personal Radiation Protection](#)

**Tag:** [Lead Gloves - Features hook and loop for full dexterity per pair.](#)



## PEG RACK WITH 8 PEGS

- Peg Rack is economical for the quick and easy storage of single or multiple aprons.
- Powder coated steel construction.
- Easy installation.

[Read More](#)

**CODE** C-PRO1008

**Category:** [Personal Radiation Protection](#)

**Tags:** [Hangers](#), [lead shielding](#), [Racks & Shielding](#)





## VINYL LEAD SEAMLESS GLOVES

[Read More](#)


**CODE** C-PRO1015

**Category:** [Personal Radiation Protection](#)

**Tag:** [Lead Gloves - Seamless lead vinyl gloves with foam liner \(available in Dark blue or Burgundy\) per pair](#)



## COLLIMATOR FILTERS



### LATERAL STERNUM

Lateral Sternum: Shield for lateral exposures

[Read More](#)

**CODE** NFSB

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



## NOLAN FILTERS, A-P CERVICO THORACIC

A-P Cervico Thoracic Compensating wedge

[Read More](#)

**CODE** NFAPCT

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#)





## NOLAN FILTER POINT FILTERS

***POINT FILTERS. Can be used for compensation in all views.***


Adjust the exposure to the place according to anatomical thickness.  
Lower until the bottom is just above the iliac crests.

[Read More](#)

**CODE** NFPP

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)



## OCTOSTOP HALF INGOT FILTER

Half Ingot Filter

L29 x W19 x H4.5 cm approx 1.8kg.

[Read More](#)

**CODE** C-FIL1034

**Category:** [Collimator Filters](#)

**Tags:** [boomerang](#), [compensating filter](#), [filter](#), [Filters and Grids](#), [ingot](#)



# PERSONAL RADIATION PROTECTION

## RADIATION PROTECTION BARRIER

Mobile or fixed lead barriers with Lead glass window - manufactured to size.

[Read More](#)

**CODE** N/A

**Categories:** [Lead Glass and Mobile Shields](#), [Personal Radiation Protection](#)






## CLOSED LOOP APRON HANGER

[Read More](#)

**CODE** C-PRO1014

**Categories:** [Personal Radiation Protection](#), [Radiation Protection](#)

**Tags:** [Apron Hanger](#), [Hanger](#), [Radiation Protection](#), [Steel Hanger](#)



## HEAVY-DUTY APRON HANGER

[Read More](#)

**CODE** C-PRO1013

**Categories:** [Personal Radiation Protection](#), [Radiation Protection](#)

**Tags:** [Hanger](#), [Lead Gown](#), [Radiation Protection](#), [Steel Hanger](#)



## WALL MOUNT SWING ARM APRON RACK

Wall Mount Swing Arm Apron Rack features 5 individual chrome shoulder shape swing arms that fold flush against the wall. Available in right and left folding models. Simply hang two models side-by-side for a 10 hanger wall console unit.

[Read More](#)

**CODE** C-PRO1005

**Category:** [Personal Radiation Protection](#)

**Tags:** [Gown](#), [Lead Gown](#), [Lead gown rack](#), [Rack](#)

## LEAD THYROID COLLAR

**Quick Release Collar** offers the same look and comfort as our model 621U, but without the hook & loop closure. Rear buckle closure. Available in our wide selection of nylons, ripstops and designer fabrics.

- **Standard Protection:** 0.5mm Pb Equivalency

**New Style Thyroid Collar** offers a trim yet shielded and comfortable design. Includes binding and rear hook & loop closure.

- **Standard Protection:** 0.5mm Pb Equivalency

Please note: Colour is light blue

[Read More](#)

**CODE** C-PRO1003-L

**Category:** [Personal Radiation Protection](#)

**Tags:** [lead protection](#), [Thyroid Collar](#)







## LEAD GOWN

Choose your Size

(XS) X-Small: 18" / 46 cm Wide - Fits Chest 30-34" / 76.2-86.4cm

(S) Small: 20" / 51 cm Wide - Fits Chest 34-38" / 86.4-96.5cm

(M) Medium: 22" / 56 cm Wide - Fits chest 38-42" / 96.5-106.7cm

(L) Large: 24" / 61 cm Wide - Fits Chest 42-46" 106.7-116.8cm


Regular lead = 0.5mm

[Read More](#)

**CODE** C-PRO1000

**Category:** [Personal Radiation Protection](#)

**Tags:** [Lead Gown](#), [lead shielding](#), [xray](#), [xray gown](#)



## MOBILE APRON RACK 1.17MT

[Read More](#)

**CODE** C-PRO1009

**Category:** [Personal Radiation Protection](#)

**Tags:** [Hangers](#), [Racks & Shielding](#)



## VET MITT

[Read More](#)

**CODE** C-PRO1002

**Category:** [Personal Radiation Protection](#)

**Tag:** [Lead Gloves](#)



## OPEN PALM SHIELD

[Read More](#)

**CODE** C-PRO1001

**Category:** [Personal Radiation Protection](#)

**Tag:** [Lead Gloves - Features hook and loop for full dexterity per pair.](#)



## PEG RACK WITH 8 PEGS

- Peg Rack is economical for the quick and easy storage of single or multiple aprons.
- Powder coated steel construction.
- Easy installation.

[Read More](#)

**CODE** C-PRO1008

**Category:** [Personal Radiation Protection](#)

**Tags:** [Hangers](#), [lead shielding](#), [Racks & Shielding](#)





## VINYL LEAD SEAMLESS GLOVES

[Read More](#)

**CODE** C-PRO1015

**Category:** [Personal Radiation Protection](#)

**Tag:** [Lead Gloves - Seamless lead vinyl gloves with foam liner \(available in Dark blue or Burgundy\) per pair](#)

