

# COLLIMATOR FILTERS

## LATERAL STERNUM

Lateral Sternum: Shield for lateral exposures

[Read More](#)

**CODE** NFSB

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Nolan Filters](#)





## NOLAN FILTERS, A-P CERVICO THORACIC

A-P Cervico Thoracic Compensating wedge

[Read More](#)

**CODE** NFAPCT

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#)



## NOLAN FILTER POINT FILTERS

***POINT FILTERS. Can be used for compensation in all views.***

Adjust the exposure to the place according to anatomical thickness.  
Lower until the bottom is just above the iliac crests.

[Read More](#)

**CODE** NFPP

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)



## OCTOSTOP HALF INGOT FILTER

Half Ingot Filter

L29 x W19 x H4.5 cm approx 1.8kg.

[Read More](#)

**CODE** C-FIL1034

**Category:** [Collimator Filters](#)

**Tags:** [boomerang](#), [compensating filter](#), [filter](#), [Filters and Grids](#), [ingot](#)

