

COLLIMATOR FILTERS



LATERAL STERNUM

Lateral Sternum: Shield for lateral exposures

[Read More](#)

CODE NFSB

Category: [Collimator Filters](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTERS, A-P CERVICO THORACIC

A-P Cervico Thoracic Compensating wedge

[Read More](#)

CODE NFAPCT

Category: [Collimator Filters](#)

Tags: [Filters](#), [Full Spine Filters](#)



NOLAN FILTER POINT FILTERS

POINT FILTERS. Can be used for compensation in all views.


Adjust the exposure to the place according to anatomical thickness.
Lower until the bottom is just above the iliac crests.

[Read More](#)

CODE NFPP

Category: [Collimator Filters](#)

Tags: [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)



OCTOSTOP HALF INGOT FILTER

Half Ingot Filter

L29 x W19 x H4.5 cm approx 1.8kg.

[Read More](#)

CODE C-FIL1034

Category: [Collimator Filters](#)

Tags: [boomerang](#), [compensating filter](#), [filter](#), [Filters and Grids](#), [ingot](#)

