

PERSONAL RADIATION PROTECTION

RADIATION PROTECTION BARRIER

Mobile or fixed lead barriers with Lead glass window - manufactured to size.

[Read More](#)

CODE N/A

Categories: [Lead Glass and Mobile Shields](#), [Personal Radiation Protection](#)






CLOSED LOOP APRON HANGER

[Read More](#)

CODE C-PRO1014

Categories: [Personal Radiation Protection](#), [Radiation Protection](#)

Tags: [Apron Hanger](#), [Hanger](#), [Radiation Protection](#), [Steel Hanger](#)



HEAVY-DUTY APRON HANGER

[Read More](#)

CODE C-PRO1013

Categories: [Personal Radiation Protection](#), [Radiation Protection](#)

Tags: [Hanger](#), [Lead Gown](#), [Radiation Protection](#), [Steel Hanger](#)



WALL MOUNT SWING ARM APRON RACK

Wall Mount Swing Arm Apron Rack features 5 individual chrome shoulder shape swing arms that fold flush against the wall. Available in right and left folding models. Simply hang two models side-by-side for a 10 hanger wall console unit.

[Read More](#)

CODE C-PRO1005

Category: [Personal Radiation Protection](#)

Tags: [Gown](#), [Lead Gown](#), [Lead gown rack](#), [Rack](#)



LEAD THYROID COLLAR

Quick Release Collar offers the same look and comfort as our model 621U, but without the hook & loop closure. Rear buckle closure. Available in our wide selection of nylons, ripstops and designer fabrics.

- **Standard Protection:** 0.5mm Pb Equivalency

New Style Thyroid Collar offers a trim yet shielded and comfortable design. Includes binding and rear hook & loop closure.

- **Standard Protection:** 0.5mm Pb Equivalency

Please note: Colour is light blue

[Read More](#)

CODE C-PRO1003-L

Category: [Personal Radiation Protection](#)

Tags: [lead protection](#), [Thyroid Collar](#)



LEAD GOWN

Choose your Size

(XS) X-Small: 18" / 46 cm Wide - Fits Chest 30-34" / 76.2-86.4cm

(S) Small: 20" / 51 cm Wide - Fits Chest 34-38" / 86.4-96.5cm

(M) Medium: 22" / 56 cm Wide - Fits chest 38-42" / 96.5-106.7cm

(L) Large: 24" / 61 cm Wide - Fits Chest 42-46" 106.7-116.8cm


Regular lead = 0.5mm

[Read More](#)

CODE C-PRO1000

Category: [Personal Radiation Protection](#)

Tags: [Lead Gown](#), [lead shielding](#), [xray](#), [xray gown](#)



MOBILE APRON RACK 1.17MT

[Read More](#)

CODE C-PRO1009

Category: [Personal Radiation Protection](#)

Tags: [Hangers](#), [Racks & Shielding](#)



VET MITT

[Read More](#)

CODE C-PRO1002

Category: [Personal Radiation Protection](#)

Tag: [Lead Gloves](#)



OPEN PALM SHIELD

[Read More](#)

CODE C-PRO1001

Category: [Personal Radiation Protection](#)

Tag: [Lead Gloves - Features hook and loop for full dexterity per pair.](#)



PEG RACK WITH 8 PEGS

- Peg Rack is economical for the quick and easy storage of single or multiple aprons.
- Powder coated steel construction.
- Easy installation.

[Read More](#)

CODE C-PRO1008

Category: [Personal Radiation Protection](#)

Tags: [Hangers](#), [lead shielding](#), [Racks & Shielding](#)





VINYL LEAD SEAMLESS GLOVES

[Read More](#)

CODE C-PRO1015

Category: [Personal Radiation Protection](#)

Tag: [Lead Gloves - Seamless lead vinyl gloves with foam liner \(available in Dark blue or Burgundy\) per pair](#)

