

# PERSONAL RADIATION PROTECTION



## RADIATION PROTECTION BARRIER

Mobile or fixed lead barriers with Lead glass window - manufactured to size.

[Read More](#)

**CODE** N/A

**Categories:** [Lead Glass and Mobile Shields](#), [Personal Radiation Protection](#)




## CLOSED LOOP APRON HANGER

[Read More](#)

**CODE** C-PRO1014

**Categories:** [Personal Radiation Protection](#), [Radiation Protection](#)

**Tags:** [Apron Hanger](#), [Hanger](#), [Radiation Protection](#), [Steel Hanger](#)



## HEAVY-DUTY APRON HANGER

[Read More](#)

**CODE** C-PRO1013

**Categories:** [Personal Radiation Protection](#), [Radiation Protection](#)

**Tags:** [Hanger](#), [Lead Gown](#), [Radiation Protection](#), [Steel Hanger](#)



## WALL MOUNT SWING ARM APRON RACK

Wall Mount Swing Arm Apron Rack features 5 individual chrome shoulder shape swing arms that fold flush against the wall. Available in right and left folding models. Simply hang two models side-by-side for a 10 hanger wall console unit.

[Read More](#)

**CODE** C-PRO1005

**Category:** [Personal Radiation Protection](#)

**Tags:** [Gown](#), [Lead Gown](#), [Lead gown rack](#), [Rack](#)

## LEAD THYROID COLLAR

**Quick Release Collar** offers the same look and comfort as our model 621U, but without the hook & loop closure. Rear buckle closure. Available in our wide selection of nylons, ripstops and designer fabrics.

- **Standard Protection:** 0.5mm Pb Equivalency

**New Style Thyroid Collar** offers a trim yet shielded and comfortable design. Includes binding and rear hook & loop closure.

- **Standard Protection:** 0.5mm Pb Equivalency

Please note: Colour is light blue

[Read More](#)

**CODE** C-PRO1003-L

**Category:** [Personal Radiation Protection](#)

**Tags:** [lead protection](#), [Thyroid Collar](#)





## LEAD GOWN

Choose your Size

(XS) X-Small: 18" / 46 cm Wide - Fits Chest 30-34" / 76.2-86.4cm

(S) Small: 20" / 51 cm Wide - Fits Chest 34-38" / 86.4-96.5cm

(M) Medium: 22" / 56 cm Wide - Fits chest 38-42" / 96.5-106.7cm

(L) Large: 24" / 61 cm Wide - Fits Chest 42-46" 106.7-116.8cm


Regular lead = 0.5mm

[Read More](#)

**CODE** C-PRO1000

**Category:** [Personal Radiation Protection](#)

**Tags:** [Lead Gown](#), [lead shielding](#), [xray](#), [xray gown](#)



## MOBILE APRON RACK 1.17MT

[Read More](#)

**CODE** C-PRO1009

**Category:** [Personal Radiation Protection](#)

**Tags:** [Hangers](#), [Racks & Shielding](#)



## VET MITT

[Read More](#)

**CODE** C-PRO1002

**Category:** [Personal Radiation Protection](#)

**Tag:** [Lead Gloves](#)





## OPEN PALM SHIELD

[Read More](#)

**CODE** C-PRO1001

**Category:** [Personal Radiation Protection](#)

**Tag:** [Lead Gloves - Features hook and loop for full dexterity per pair.](#)



## PEG RACK WITH 8 PEGS

- Peg Rack is economical for the quick and easy storage of single or multiple aprons.
- Powder coated steel construction.
- Easy installation.

[Read More](#)

**CODE** C-PRO1008

**Category:** [Personal Radiation Protection](#)

**Tags:** [Hangers](#), [lead shielding](#), [Racks & Shielding](#)





## VINYL LEAD SEAMLESS GLOVES

[Read More](#)

**CODE** C-PRO1015

**Category:** [Personal Radiation Protection](#)

**Tag:** [Lead Gloves - Seamless lead vinyl gloves with foam liner \(available in Dark blue or Burgundy\) per pair](#)

