

MARKERS

BEEKLEY CT-SPOT® 123

CT-SPOT® 123

Ref: 123, (QTY 50/box)

Ideal for marking 3-point set-ups, isocenters, and field borders, these unique markers image brightly during simulation for visibility through multiple slices.


- 2.0cm x 2.0cm non-metallic crosshair
- for 3-point set-ups, isocenters, and field borders

[Read More](#)

CODE BEEK123

Category: [Markers](#)

Tags: [CT Simulation Skin Markers](#), [Grids and Markers](#)



BEEKLEY CT-SPOT® 120

Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.

BEEKLEY CT-SPOT® 120, Ref: 120 (QTY: 85/box)


- 4.0mm non-metallic pellet
- for 3-point set-ups and isocenters

[Read More](#)

CODE C-BEEK120

Categories: [Markers](#), [Skin Markers](#)

Tag: [CT Simulation Skin Markers](#)



BEEKLEY CT-SPOT® 119

BEEKLEY CT-SPOT® 119

Ref: 119, (QTY: 85/box)

Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.


- 2.3mm non-metallic pellet
- for 3-points set-ups and isocenters

[Read More](#)

CODE C-BEEK119

Categories: [Markers](#), [Skin Markers](#)

Tags: [CT Simulation Skin Markers](#), [Grids and Markers](#)





BEEKLEY CT-SPOT® 118

CT-SPOT® 118

**Ref: 118, (QTY: 320cm cut to measure
roll/box)**

- 2.0mm diameter non-metallic line
- for marking scars and field borders

[Read More](#)

CODE BEEK118

Category: [Markers](#)

Tag: [CT Simulation Skin Markers](#)



SKIN MARKERS

BEEKLEY CT-SPOT® 120

Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.

BEEKLEY CT-SPOT® 120, Ref: 120 (QTY: 85/box)


- 4.0mm non-metallic pellet
- for 3-point set-ups and isocenters

[Read More](#)

CODE C-BEEK120

Categories: [Markers](#), [Skin Markers](#)

Tag: [CT Simulation Skin Markers](#)



BEEKLEY CT-SPOT® 119

BEEKLEY CT-SPOT® 119

Ref: 119, (QTY: 85/box)

Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.

- 2.3mm non-metallic pellet
- for 3-points set-ups and isocenters

[Read More](#)

CODE C-BEEK119

Categories: [Markers](#), [Skin Markers](#)

Tags: [CT Simulation Skin Markers](#), [Grids and Markers](#)

