

POSITIONING AIDS



POSITIONING FOAM, FOREHEAD IMMOBILISER 120X270X250MM

Forehead Immobiliser:

Is a U shaped designed to support the forehead of a patient during X-ray, CT and MRI studies. Can be covered in Dartex or vinyl. 120x270x250mm.

[Read More](#)

CODE C-FOAM62

Category: [Foams](#)

Tag: [Forehead Immobiliser](#)



FOAMS

Imaging Positioning, Patient Positioning, Support Foams,

Manufactured by ATX



POSITIONING FOAM, FOREHEAD IMMOBILISER 120X270X250MM

Forehead Immobiliser:

Is a U shaped designed to support the forehead of a patient during X-ray, CT and MRI studies. Can be covered in Dartex or vinyl. 120x270x250mm.

[Read More](#)

CODE C-FOAM62

Category: [Foams](#)

Tag: [Forehead Immobiliser](#)



