

# COLLIMATOR FILTERS

## FULL SET SPINE FILTERS


Complete box set of all Full Spine Filters. **Reduce Radiation dose by up to 95%.**

[Read More](#)

**CODE** C-FIL1025


**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Filters and Grids](#), [Full Spine Filters](#)





# RADIATION PROTECTION



## NOLAN FILTER FACE PLATE

Mounts directly on to the collimator to  
accommodate multiple filters

[Read More](#)

**CODE** NFFP

**Category:** [Radiation Protection](#)

**Tags:** [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)



## NOLAN FILTERS, A-P CERVICO THORACIC

A-P Cervico Thoracic Compensating wedge

[Read More](#)

**CODE** NFAPCT

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#)



## NOLAN FILTER POINT FILTERS

***POINT FILTERS. Can be used for compensation in all views.***

Adjust the exposure to the place according to anatomical thickness.  
Lower until the bottom is just above the iliac crests.

[Read More](#)

**CODE** NFPP

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)



## COLLIMATOR FILTERS



### NOLAN FILTERS, A-P CERVICO THORACIC

A-P Cervico Thoracic Compensating wedge

[Read More](#)

**CODE** NFAPCT

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#)



## NOLAN FILTER POINT FILTERS

***POINT FILTERS. Can be used for compensation in all views.***

Adjust the exposure to the place according to anatomical thickness.  
Lower until the bottom is just above the iliac crests.

[Read More](#)

**CODE** NFPP

**Category:** [Collimator Filters](#)

**Tags:** [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)





