

COLLIMATOR FILTERS

FULL SET SPINE FILTERS


Complete box set of all Full Spine Filters. **Reduce Radiation dose by up to 95%.**

[Read More](#)


CODE C-FIL1025

Category: [Collimator Filters](#)

Tags: [Filters](#), [Filters and Grids](#), [Full Spine Filters](#)



RADIATION PROTECTION



NOLAN FILTER FACE PLATE

Mounts directly on to the collimator to
accommodate multiple filters

[Read More](#)

CODE NFFP

Category: [Radiation Protection](#)

Tags: [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)



NOLAN FILTERS, A-P CERVICO THORACIC

A-P Cervico Thoracic Compensating wedge

[Read More](#)

CODE NFAPCT

Category: [Collimator Filters](#)

Tags: [Filters](#), [Full Spine Filters](#)



NOLAN FILTER POINT FILTERS

POINT FILTERS. Can be used for compensation in all views.

Adjust the exposure to the place according to anatomical thickness.
Lower until the bottom is just above the iliac crests.

[Read More](#)

CODE NFPP

Category: [Collimator Filters](#)

Tags: [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)

COLLIMATOR FILTERS



NOLAN FILTERS, A-P CERVICO THORACIC

A-P Cervico Thoracic Compensating wedge

[Read More](#)

CODE NFAPCT

Category: [Collimator Filters](#)

Tags: [Filters](#), [Full Spine Filters](#)



NOLAN FILTER POINT FILTERS

POINT FILTERS. Can be used for compensation in all views.

Adjust the exposure to the place according to anatomical thickness.
Lower until the bottom is just above the iliac crests.

[Read More](#)

CODE NFPP

Category: [Collimator Filters](#)

Tags: [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)

