

MARKERS

BEEKLEY CT-SPOT® 123

CT-SPOT® 123

Ref: 123, (QTY 50/box)

Ideal for marking 3-point set-ups, isocenters, and field borders, these unique markers image brightly during simulation for visibility through multiple slices.


- 2.0cm x 2.0cm non-metallic crosshair
- for 3-point set-ups, isocenters, and field borders

[Read More](#)

CODE BEEK123

Category: [Markers](#)

Tags: [CT Simulation Skin Markers](#), [Grids and Markers](#)






BEEKLEY CT-SPOT® 119

BEEKLEY CT-SPOT® 119

Ref: 119, (QTY: 85/box)

Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.

- 2.3mm non-metallic pellet
- for 3-points set-ups and isocenters



[Read More](#)

CODE C-BEEK119

Categories: [Markers](#), [Skin Markers](#)

Tags: [CT Simulation Skin Markers](#), [Grids and Markers](#)

SKIN MARKERS

BEEKLEY CT-SPOT® 119

BEEKLEY CT-SPOT® 119

Ref: 119, (QTY: 85/box)

Ideal for marking 3-point set-ups and isocenters, these unique markers image brightly during simulation without streaking or artifact.

- 2.3mm non-metallic pellet
- for 3-points set-ups and isocenters

[Read More](#)

CODE C-BEEK119

Categories: [Markers](#), [Skin Markers](#)

Tags: [CT Simulation Skin Markers](#), [Grids and Markers](#)

