

MARKERS

BEEKLEY MARK-THROUGH GUIDELINES® 219

Mark-Through GuideLines® 219

Ref: 219, (QTY: 15/box)

- 6.25" (15.88cm) x 6.75" (17.15cm) biopsy grid
- 11 non-lead radiopaque lines
- mark-through material for easy mark placement anywhere on grid
- ideal for CT-guided biopsies, drainage procedures, and pain management injections

[Read More](#)

CODE BEEK219

Category: [Markers](#)

Tag: [Mark-Through GuideLines](#)



