

RADIATION PROTECTION



NOLAN FILTER, GONAD SHIELD SYSTEM SET

Set of male and female gonad shields
and belt.

[Read More](#)

CODE NFAPSET

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER, GONAD SYSTEM BELT

Velcro belt to suspend gonad shields

[Read More](#)

CODE NFAPSBELT

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER, A-P FEMALE GONAD SHIELD

A-P Female Gonad Shield for AP shots.

[Read More](#)

CODE NFAPFGS

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER, A-P MALE GONAD

A-P Male Gonad Shield for AP shots.

[Read More](#)

CODE NFAPMGS

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER, LATERAL SPINOUS

Lateral Spinous: Reduces scatter and skin dose and
balances plate dose by blocking x-rays to the spinous process.

[Read More](#)

CODE NFLS

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER, A-P PARASPINAL BLOCKER (TOTAL)

A-P Paraspinal Blocker (Total): Totally blocks anatomy either side of the spine.

[Read More](#)

CODE NFAPSBF

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER, A-P PARASPINAL BLOCKER (PARTIAL)

A-P Paraspinal Blocker (Partial): Partially blocks anatomy either side of the spine.

[Read More](#)

CODE NFAPSBP

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER, LATERAL GONAD

Lateral M/F Gonad shield for lateral exposures

[Read More](#)

CODE NFLG

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER, LATERAL LUMBAR LUNG

Lateral Lumbar Lung: Lumbar Lung shield for lateral exposures

[Read More](#)

CODE NFLLL

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



LATERAL STERNUM

Lateral Sternum: Shield for lateral exposures

[Read More](#)

CODE NFSB

Category: [Collimator Filters](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER, LATERAL THORACIC

Lateral Thoracic. Compensating wedge

[Read More](#)

CODE NFLT

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER LATERAL CERVICO THORACIC


Accommodates relative thickness of the shoulder to the Cervical or Thyroid region.

[Read More](#)

CODE NFLCT

Category: [Radiation Protection](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER FACE PLATE

Mounts directly on to the collimator to
accommodate multiple filters

[Read More](#)

CODE NFFP

Category: [Radiation Protection](#)

Tags: [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)



NOLAN FILTER POINT FILTERS

POINT FILTERS. Can be used for compensation in all views.

Adjust the exposure to the place according to anatomical thickness.
Lower until the bottom is just above the iliac crests.

[Read More](#)

CODE NFPP

Category: [Collimator Filters](#)

Tags: [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)

COLLIMATOR FILTERS



LATERAL STERNUM

Lateral Sternum: Shield for lateral exposures

[Read More](#)

CODE NFSB

Category: [Collimator Filters](#)

Tags: [Filters](#), [Filters and Grids](#), [Nolan Filters](#)



NOLAN FILTER POINT FILTERS

POINT FILTERS. Can be used for compensation in all views.

Adjust the exposure to the place according to anatomical thickness.
Lower until the bottom is just above the iliac crests.

[Read More](#)

CODE NFPP

Category: [Collimator Filters](#)

Tags: [Filters](#), [Full Spine Filters](#), [Nolan Filters](#)

