

Amplitude	1 mV
Acuity Tube Current	2 mA
Focal Spot Size	0.4 x 0.4 mm
Output Power	142 mW
Exposure Time	0.01 - 1.3 sec

Mechanical

Anode Angle	12.5 degrees
Inherent Filtration	At least 1.5 mmAl
Anode Coverage	±160 mm at 180-250 mm

Dimension

	165x2 x 105-266x5 x 153x6 (mm/inch)
--	-------------------------------------

Weight

	About 1.7 kg
--	--------------

Power requirement

Rechargeable battery

REMEX-T100 PORTABLE X-RAY DIGITAL CAMERA

***Focal Spot Value : 0.4x0.4mm.**

Tags: [dental](#), [Remedi](#), [Remex](#), [T100](#)



