

X-RAY EQUIPMENT

REMEX-T100 PORTABLE X-RAY DIGITAL CAMERA

Remex-T100 handheld portable X-ray generator

- *Tube Voltage : 70Kv
- *Tube Current : 2mA
- *Focal Spot Value : 0.4x0.4mm.

[Read More](#)
CODE E-GEN1058
Category: [X-Ray Components](#)
Tags: [dental](#), [Remedi](#), [Remex](#), [T100](#)



X-RAY COMPONENTS

REMEX-T100 PORTABLE X-RAY DIGITAL CAMERA

Remex-T100 handheld portable X-ray generator

***Tube Voltage : 70Kv**

***Tube Current : 2mA**

***Focal Spot Value : 0.4x0.4mm.**

[Read More](#)
CODE E-GEN1058
Category: [X-Ray Components](#)
Tags: [dental](#), [Remedi](#), [Remex](#), [T100](#)



Tube voltage range	20-70KV
Tube Voltage Control	0-70KV
Focal Spot Size	0.4x0.4mm
Exposure Time	10-99sec
Exposure Time	0.01 - 1.5min

Mechanical	
Tube Angle	15.5 degrees
Inertial Motor	About 1.5 sec/0.5
Auto Storage	100 sec at 80-100 sec

Dimensions	190x110x200x110x110mm (30x40x110x40x110mm)
Weight	About 1.7kg

Power requirement	
Rechargeable battery	

