

1800 330 118

ATX

VETERINARY SOLUTIONS

www.atxvet.com.au

ATX PORTABLE X-RAY

This exclusive to ATX system re-defines portable X-ray systems. With a smart interface options enables an unprecedented amount of information at your finger tips.



	ATX HF16	ATX SHF24	ATX SHF32
Max Output Power	1.6KW	2.4KW	3.2KW
Generator Frequency	100KHz	70KHz	70KHz
kV Range	40~90KV	40~100KV	40~100KV
mA Range	12~30mA	16~40mA	25~60mA
mAs Range	0.4~100mAs	0.4~100mAs	0.4~100mAs
Focal Spot	1.2mm * 1.2mm	1.2mm * 1.2mm	2.0mm * 2.0mm
Target Angle	19 °	16 °	16 °
Anode Heat Storage	10KHU	50KHU	40KHU
Input Voltage	110V / 220V	110V / 220V	200V / 240V
Phase/Frequency	Single / 50~60Hz	Single / 50~60Hz	Single / 50~60Hz
Dimension(mm)	335*153*330	360*200*190	360*200*190
Packing Size(mm)	520*400*330	520*400*330	520*400*330
Net weight	6.8KG	11KG	12KG
G-weight	20KG	22KG	22KG
Smart Interface		yes	yes

ATX Smart Portable X-ray machine delivers many benefits and features. Specific and comprehensive APR's (Anatomical Presets) can be loaded on each machine to suit the practice and acquisition system, giving you unique control. A log of the number of x-rays and settings of each exposure is recorded, so this can be matched against invoiced studies. Remote access is possible for fast diagnosis, giving you access when you need it most and also lets us support you with fast assistance. Built in calibration and performance tools assist to verify the operation of the machine. This simple and efficient machine will not disappoint.

www.atxvet.com.au

1800 330 118

vet@atx.com.au